

PRECISION - STRENGTH - HEALING



SanaCuraTM
SUTURES

By Axa Parenterals Ltd.

We Comply with all major certifications,
Including ISO 9001:2015, EN ISO 13485:2016, CE, and GMP

GLOBAL PARTNERSHIP FOR SURGICAL EXCELLENCE

Notes:





About Us:

SanaCura Sutures is a specialized vertical under the trusted leadership of Axa Parenterals Ltd. and Heilsa Life Sciences Pvt. Ltd., dedicated to delivering high-quality surgical sutures for healthcare professionals worldwide.

Rooted in excellence, SanaCura upholds the same commitment to innovation, precision, and quality that has defined our parent companies for decades. With strict quality control systems in place, we have quickly gained recognition as a reliable manufacturer, supplier, and exporter of sutures in India and across the globe.

Our state-of-the-art manufacturing facility in Uttarakhand is fully equipped with the latest machinery and technology, ensuring seamless operations and consistent product excellence.

We are driven by a clear mission—to become the most trusted suture manufacturer in India—and a strong vision—to make surgical solutions more affordable and accessible without compromising on quality.

Backed by:

- Axa Parenterals Ltd., a leader in sterile parenteral preparations serving over 18,000 hospitals across 25+ states with more than 100 product launches.

- Heilsa Life Sciences Pvt. Ltd., a company known for its EU-GMP compliant manufacturing of critical care injectables including liquid vials, ampoules, pre-filled syringes, and more.

Together, we bring deep pharmaceutical expertise, robust infrastructure, and a passion for improving healthcare outcomes through reliable surgical solutions.



SanaSynth

Braided & Coated Polyglactin 910
Absorbable Surgical Suture U.S.P.
Sterile, Synthetic Surgical Suture



SanaSynth is a synthetic, absorbable, sterile suture made from 90% glycolide and 10% L-lactide, ideal for soft tissue closure, especially skin and mucosa.

It offers strong knot security, smooth texture for minimal trauma, and superior tensile strength supporting wounds for 28–32 days with full absorption in 56–70 days.

Applications include uterine closure and soft tissue approximation. It is a braided, synthetic absorbable suture.

Specifications:

| | |
|--------------------|--|
| Type | : Braided, Synthetic Absorbable |
| Composition | : Polyglactin 910 |
| Coating | : Poly (glycolide-co-L-lactide) & Calcium Stearate |
| Color | : Violet |
| Strength Retention | : 75% (14 days), 50% (21 days) |
| Absorption | : 56–70 days via hydrolysis |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 5 years |
| Packaging | : Aluminum foil packs (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

SanaSynth Fast

Braided & Coated Polyglactin 910 - Fast
Absorbable Surgical Suture U.S.P.
Sterile, Synthetic Surgical Suture



SanaSynth Fast is a synthetic, absorbable, sterile suture made from 90% glycolide and 10% L-lactide, designed for soft tissue closure with short-term wound support. It offers strong knot security, smooth texture, and superior tensile strength that supports wounds for 28–32 days, with complete absorption in 28–42 days.

Applications include uterine closure and soft tissue approximation. It is a braided, synthetic absorbable suture.

Specifications:

| | |
|--------------------|--|
| Type | : Braided, Synthetic Absorbable |
| Composition | : Polyglactin 910 |
| Coating | : Poly (glycolide-co-L-lactide) & Calcium Stearate |
| Color | : Undyed |
| Strength Retention | : 75% (3 days), 50% (5 days) |
| Absorption | : 28–42 days via hydrolysis |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 3 years |
| Packaging | : Aluminum foil packs (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

| Product | Product Code | Description | Units / Box |
|-----------------------|----------------|--|-------------|
| Polyglactin 910 | MS 2478 | USP 2 90 cm 1/2 Circle Reverse Cutting (H) 40 mm | 12 Foils |
| Polyglactin 910 | MS 2347 | USP 1 90 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglactin 910 | MS 2347VS | USP 1 110 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglactin 910 | MS 2347N | USP 1 115 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglactin 910 | MS 2347VSS | USP 1 125 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglactin 910 | MS 2347DA | USP 1 140 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 | MS 2347DA/180 | USP 1 180 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 | MS 2347DA/180N | USP 1 180 cm 1/2 Circle Round Body 50 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 | MS 2421 | USP 1 90 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglactin 910 | MS 2448 | USP 1 100 cm 1/2 Circle Round Body (H) 50 mm | 12 Foils |
| Polyglactin 910 | MS 2350RB | USP 1 120 cm 1/2 Circle Round Body 36 mm | 12 Foils |
| Polyglactin 910 | MS 2360 | USP 1 90 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglactin 910 | MS 2363 | USP 1 70 cm 1/2 Circle Taper Cutting (H) 40 mm | 12 Foils |
| Polyglactin 910 | MS 2438N | USP 1 90 cm 1/2 Circle Reverse Cutting 45 mm | 12 Foils |
| Polyglactin 910 | MS 2519 | USP 1 90 cm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 | MS 2826 | USP 1 70 cm 1/2 Circle Reverse Cutting 23 mm | 12 Foils |
| Polyglactin 910 | MS 2826A | USP 1 35 cm 1/2 Circle Reverse Cutting 23 mm | 12 Foils |
| Polyglactin 910 | MS 2438 | USP 1 90 cm 1/2 Circle Reverse Cutting (H) 50 mm | 12 Foils |
| Polyglactin 910 | MS 2348 | USP 1 75 cm 1/2 Circle Round Body 48 mm | 12 Foils |
| Polyglactin 910 - (U) | MS 947 | USP 1 90 cm 1/2 Circle Round Body 36 mm | 12 Foils |
| Polyglactin 910 | MS 2346 | USP 0 90 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglactin 910 | MS 2346VS | USP 0 110 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglactin 910 | MS 2346DA | USP 0 140 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 | MS 2346DA/180 | USP 0 180 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 | MS 2534 | USP 0 90 cm 1/2 Circle Reverse Cutting 36 mm | 12 Foils |
| Polyglactin 910 | MS 2518 | USP 0 90 cm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 | MS 2517OB | USP 0 90 cm 1/2 Circle Taper Cutting 36 mm | 12 Foils |
| Polyglactin 910 | MS 2518DA/180 | USP 0 180 cm 1/2 Circle Taper Cutting 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 - (U) | MS 946 | USP 0 90 cm 1/2 Circle Round Body 36 mm | 12 Foils |
| Polyglactin 910 | MS 2338 | USP 0 90 cm 1/2 Circle Round Body 30 mm | 12 Foils |

| Product | Product Code | Description | Units / Box |
|-----------------------|---------------|---|-------------|
| Polyglactin 910 | MS 2825 | USP 0 45 cm 1/2 Circle Reverse Cutting 23 mm | 12 Foils |
| Polyglactin 910 | MS 2317 | USP 2-0 90 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglactin 910 | MS 2317DA | USP 2-0 140 cm 1/2 Circle Round Body 30 mm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polyglactin 910 | MS 2345 | USP 2-0 90 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglactin 910 | MS 2382 | USP 2-0 90 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglactin 910 | MS 2404 | USP 2-0 90 cm 3/8 Circle Conventional Cutting 26 mm | 12 Foils |
| Polyglactin 910 - (U) | MS 945 | USP 2-0 90 cm 1/2 Circle Round Body 36 mm | 12 Foils |
| Polyglactin 910 | MS 2341 | USP 2-0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglactin 910 - (U) | MS 2404 U | USP 2-0 90 cm 3/8 Circle Cutting 40 mm | 12 Foils |
| Polyglactin 910 | MS 9828 | USP 2-0 75 cm 1/2 Circle Taper cutting 22 mm | 12 Foils |
| Polyglactin 910 | MS 2345DA/180 | USP 2-0 180 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 | MS 2437 | USP 3-0 70 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Polyglactin 910 | MS 2437A | USP 3-0 70 cm 1/2 Circle Round Body 16 mm | 12 Foils |
| Polyglactin 910 | MS 2340 | USP 3-0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglactin 910 | MS 2471 | USP 3-0 45 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Polyglactin 910 | MS 2516 | USP 3-0 90 cm 1/2 Circle Taper Cutting 36 mm | 12 Foils |
| Polyglactin 910 | MS 2401 | USP 3-0 45 cm 3/8 Circle Cutting 22 mm | 12 Foils |
| Polyglactin 910 | MS 2472 | USP 3-0 90 cm 1/2 Circle Cutting 22 mm | 12 Foils |
| Polyglactin 910 - (U) | MS 2328 | USP 3-0 90 cm 3/8 Circle Precision Cutting 25 mm | 12 Foils |
| Polyglactin 910 - (U) | MS 2407 | USP 3-0 70 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Polyglactin 910 | MS 2402PS | USP 3-0 90 cm 3/8 Circle Cutting MicroPrime 19 mm | 12 Foils |
| Polyglactin 910 | MS 2304 | USP 4-0 70 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Polyglactin 910 | MS 2944 | USP 4-0 70 cm 3/8 Circle Reverse Cutting 16 mm | 12 Foils |
| Polyglactin 910 | MS 2494 | USP 4-0 45 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Polyglactin 910 | MS 2465 | USP 4-0 45 cm 5/8 Circle Cutting 16 mm | 12 Foils |
| Polyglactin 910 | MS 2305 | USP 4-0 70 cm 1/2 Circle Round Body 16 mm | 12 Foils |
| Polyglactin 910 | MS 2493PS | USP 5-0 45 cm 3/8 Circle Reverse Cutting MicroPrime 12 mm | 12 Foils |
| Polyglactin 910 | MS 2601 | USP 0 150 cm | 12 Foils |
| Polyglactin 910 | MS 2603 | USP 3-0 150 cm | 12 Foils |
| Polyglactin 910 | MS 2604 | USP 4-0 150 cm | 12 Foils |

| Product | Product Code | Description | Units / Box |
|-----------------|--------------|---|-------------|
| Polyglactin 910 | MS 2763 F | USP 0 90 cm 1/2 Circle Taper Cutting 35 mm | 12 Foils |
| Polyglactin 910 | MS 2763L F | USP 0 100 cm 1/2 Circle Taper Cutting Epicio 35 mm | 12 Foils |
| Polyglactin 910 | MS 2764 F | USP 0 75 cm 1/2 Circle Round Body 48 mm | 12 Foils |
| Polyglactin 910 | MS 2763 FRB | USP 0 100 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglactin 910 | MS 2762DY F | USP 2-0 90 cm 1/2 Circle Cutting 30 mm 1/2 Circle Cutting 30 mm | 12 Foils |
| Polyglactin 910 | MS 2762 FVS | USP 2-0 110 cm 1/2 Circle Round Body 35 mm | 12 Foils |
| Polyglactin 910 | MS 2777 F | USP 2-0 140 cm 1/2 Circle Reverse Cutting 36 mm 1/2 Circle Round Body 36 mm | 12 Foils |
| Polyglactin 910 | MS 2760 F | USP 2-0 75 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Polyglactin 910 | MS 2762 F | USP 2-0 90 cm 1/2 Circle Taper Cutting 35 mm | 12 Foils |
| Polyglactin 910 | MS 2762L F | USP 2-0 100 cm 1/2 Circle Taper Cutting Epicio 35 mm | 12 Foils |
| Polyglactin 910 | MS 2762 FRB | USP 2-0 100 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglactin 910 | MS 2761 F | USP 2-0 90 cm 1/2 Circle Round Body 38 mm | 12 Foils |
| Polyglactin 910 | MS 2739 F | USP 3-0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglactin 910 | MS 2516RB F | USP 3-0 90 cm 1/2 Circle Round Body 35 mm | 12 Foils |
| Polyglactin 910 | MS 2732 F | USP 3-0 75 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Polyglactin 910 | MS 2718 F | USP 4-0 45 cm 3/8 Circle Reverse Cutting Circum 16 mm | 12 Foils |
| Polyglactin 910 | MS 2718RB F | USP 4-0 45 cm 3/8 Circle Round Body Circum 16 mm | 12 Foils |

SanaSynth Plus Antibacterial

Braided & Coated Polyglactin 910 - Fast
Absorbable Surgical Suture U.S.P.
Sterile, Synthetic Surgical Suture



SanaSynth Fast is a synthetic, absorbable, sterile suture made from 90% glycolide and 10% L-lactide, designed for soft tissue closure with short-term wound support. It offers strong knot security, smooth texture, and superior tensile strength that supports wounds for 28–32 days, with complete absorption in 28–42 days.

Applications include uterine closure and soft tissue approximation. It is a braided, synthetic absorbable suture.

Specifications:

| | |
|--------------------|--|
| Type | : Braided, Synthetic Absorbable |
| Composition | : Polyglactin 910 |
| Coating | : Poly (glycolide-co-L-lactide) & Calcium Stearate |
| Color | : Undyed |
| Strength Retention | : 75% (3 days), 50% (5 days) |
| Absorption | : 28–42 days via hydrolysis |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 3 years |
| Packaging | : Aluminum foil packs (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

| Product | Product Code | Description | Units / Box |
|-----------------------------------|---------------|---|-------------|
| Polyglactin 910 Antibacterial | MP 2478 | USP 2 90 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2347 | USP 1 90 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2347VS | USP 1 110 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2347DA | USP 1 140 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2347DA/180 | USP 1 180 cm 1/2 Circle Taper Cutting 40 mm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2350RB L | USP 1 120 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2421 | USP 1 90 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2826 | USP 1 70 cm 1/2 Circle Reverse Cutting 23 mm | 12 Foils |
| Polyglactin 910 Antibacterial (U) | MP 947 | USP 1 90 cm 1/2 Circle Round Body 36 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2364 | USP 1 90 cm 1/2 Circle Taper Cutting (H) 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2346 | USP 0 90 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2346DA/180 | USP 0 180 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2518 | USP 0 90 cm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2534 | USP 0 90 cm 1/2 Circle Reverse Cutting 36 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2317 | USP 2-0 90 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2382 | USP 2-0 90 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglactin 910 Antibacterial (U) | MP 945 | USP 2-0 90 cm 1/2 Circle Round Body 36 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2437 | USP 3-0 70 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2472 | USP 3-0 90 cm 1/2 Circle Cutting 22 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2516 | USP 3-0 90 cm 1/2 Circle Taper Cutting 36 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2516H | USP 3-0 90 cm 1/2 Circle Taper Cutting (H) 36 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2936PS | USP 3-0 70 cm 3/8 Circle Reverse Cutting 24 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2328 | USP 3-0 90 cm 3/8 Circle Cutting 25 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2304 | USP 4-0 70 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Polyglactin 910 Antibacterial | MP 2338 | USP 0 90 cm 1/2 Circle Round Body 30 mm | 12 Foils |

SanaCryl

Braided & Coated Polyglycolic Acid
Absorbable Surgical Suture U.S.P.
Sterile, Synthetic Surgical Suture



SanaCryl is a synthetic, absorbable, sterile suture made from 100% glycolide homopolymer, designed for general soft tissue approximation and ligation. It is also suitable for ophthalmic procedures but not recommended for cardiovascular or neurological tissues. It offers strong knot security, a smooth surface for minimized tissue trauma, and superior tensile strength that supports wounds for 28–30 days, with full degradation in 70–90 days.

Applications include general soft tissue approximation and ophthalmic procedures. It is a braided, synthetic absorbable suture.

Specifications:

| | |
|--------------------|--|
| Type | : Braided, Synthetic Absorbable |
| Composition | : 100% Polyglycolic Acid (PGA) |
| Coating | : Polycaprolactone & Calcium Stearate |
| Color | : Violet (also available undyed) |
| Strength Retention | : 75% (14 days), 50% (21 days) |
| Absorption | : 70–90 days via hydrolysis |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 5 years |
| Packaging | : Aluminum foil packs (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

SanaCryl Fast

Braided & Coated Polyglycolic Acid - Fast
Absorbable Surgical Suture U.S.P.
Sterile, Synthetic Surgical Suture



SanaCryl Fast is a synthetic, absorbable, sterile suture made from 100% glycolide homopolymer, designed for soft tissue approximation and ligation. It is ideal for ophthalmic procedures but not recommended for cardiovascular or neurological tissues. It offers strong knot security, a smooth surface for minimal tissue trauma, and superior tensile strength that supports wounds for 7–10 days, with complete absorption in 40–42 days.

Applications include general soft tissue approximation and ophthalmic procedures. It is a braided, synthetic absorbable suture.

Specifications:

| | |
|--------------------|--|
| Type | : Braided, Synthetic Absorbable |
| Composition | : 100% Polyglycolic Acid (PGA) |
| Coating | : Polycaprolactone & Calcium Stearate |
| Color | : Undyed |
| Strength Retention | : 75% (3 days), 50% (5 days) |
| Absorption | : 40–42 days via hydrolysis |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 3 years |
| Packaging | : Aluminum foil packs (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

| Product | Product Code | Description | Units / Box |
|-------------------|----------------|---|-------------|
| Polyglycolic Acid | N 2228 | USP 2 90 cm 1/2 Circle Round Body 45 mm | 12 Foils |
| Polyglycolic Acid | N 2228VS | USP 2 110 cm 1/2 Circle Round Body 45 mm | 12 Foils |
| Polyglycolic Acid | N 2228DA | USP 2 140 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2478 | USP 2 90 cm 1/2 Circle Reverse Cutting (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2353 | USP 2 140 cm 1/2 Circle Reverse Cutting (H) 50 mm X 2 Nos | 12 Foils |
| Polyglycolic Acid | N 2319 | USP 1 90 cm 1/2 Circle Taper Cutting OB 36 mm | 12 Foils |
| Polyglycolic Acid | N 2343 | USP 1 20 cm 1/2 Circle Round Body (H) 30 mm | 12 Foils |
| Polyglycolic Acid | N 2347 | USP 1 90 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2347A | USP 1 45 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2347AA | USP 1 20 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2347VS | USP 1 110 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2347DA | USP 1 140 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2347DA/160 | USP 1 160 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2347DA/180 | USP 1 180 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2347DA/180RB | USP 1 180 cm 1/2 Circle Round Body 40 mm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglycolic Acid | N 2347DA/180N | USP 1 180 cm 1/2 Circle Round Body 50 mm 1/2 Circle Round Body 50 mm | 12 Foils |
| Polyglycolic Acid | N 2350 | USP 1 90 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2359 | USP 1 70 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2362 | USP 1 45 cm 1/2 Circle Taper Cutting (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2363 | USP 1 70 cm 1/2 Circle Taper Cutting (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2364 | USP 1 90 cm 1/2 Circle Taper Cutting (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2364VS | USP 1 110 cm 1/2 Circle Taper Cutting (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2421DAOR | USP 1 140 cm 1/2 Circle Reverse Cutting 40 mm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglycolic Acid | N 2438 | USP 1 90 cm 1/2 Circle Reverse Cutting (H) 50 mm | 12 Foils |
| Polyglycolic Acid | N 2448 | USP 1 100 cm 1/2 Circle Round Body (H) 50 mm | 12 Foils |
| Polyglycolic Acid | N 2513 | USP 1 90 cm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polyglycolic Acid | N 2513A | USP 1 45 cm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polyglycolic Acid | N 2513AA | USP 1 20 cm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polyglycolic Acid | N 2826 | USP 1 70 cm 1/2 Circle Reverse Cutting 23 mm | 12 Foils |
| Polyglycolic Acid | N 2826A | USP 1 35 cm 1/2 Circle Reverse Cutting 23 mm | 12 Foils |

| Product | Product Code | Description | Units / Box |
|-------------------------|--------------|---|-------------|
| Polyglycolic Acid | N 2421 | USP 1 90 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2421A | USP 1 45 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglycolic Acid - (U) | N 947 | USP 1 90 cm 1/2 Circle Round Body 36 mm | 12 Foils |
| Polyglycolic Acid | N 2347N | USP 1 115 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2421VS | USP 1 110 cm 1/2 Circle Reverse Cuttin OS Needle 40 mm | 12 Foils |
| Polyglycolic Acid | N 2338 | USP 0 90 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglycolic Acid | N 2342 | USP 0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglycolic Acid | N 2344 | USP 0 20 cm 1/2 Circle Round Body (H) 30 mm | 12 Foils |
| Polyglycolic Acid | N 2346 | USP 0 90 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglycolic Acid | N 2346A | USP 0 45 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglycolic Acid | N 2346AA | USP 0 20 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglycolic Acid | N 2346VS | USP 0 110 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglycolic Acid | N 2346DA | USP 0 140 cm 1/2 Circle Round Body 40 mm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2422 | USP 0 90 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2512 | USP 0 90 cm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polyglycolic Acid | N 2512A | USP 0 45 cm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polyglycolic Acid | N 2512AA | USP 0 20 cm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polyglycolic Acid | N 2518 | USP 0 90 cm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2518VS | USP 0 110 cm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2534 | USP 0 90 cm 1/2 Circle Reverse Cutting 36 mm | 12 Foils |
| Polyglycolic Acid | N 2546 | USP 0 70 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglycolic Acid | N 2825A | USP 0 35 cm 1/2 Circle Reverse Cutting 23 mm | 12 Foils |
| Polyglycolic Acid - (U) | N 946 | USP 0 90 cm 1/2 Circle Round Body 36 mm | 12 Foils |
| Polyglycolic Acid | N 2361 | USP 0 45 cm 1/2 Circle Taper Cutting (H) 40 mm | 12 Foils |
| Polyglycolic Acid | N 2317 | USP 2-0 90 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglycolic Acid | N 2317DA | USP 2-0 140 cm 1/2 Circle Round Body 30 mm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polyglycolic Acid | N 2341 | USP 2-0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglycolic Acid | N 2345 | USP 2-0 90 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglycolic Acid | N 2382 | USP 2-0 90 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2473 | USP 2-0 70 cm 5/8 Circle Round Body 27 mm | 12 Foils |

| Product | Product Code | Description | Units / Box |
|-------------------------|--------------|---|-------------|
| Polyglycolic Acid | N 2421 | USP 1 90 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2511 | USP 2-0 90 cm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polyglycolic Acid | N 2511A | USP 2-0 45 cm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polyglycolic Acid | N 2511AA | USP 2-0 20 cm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polyglycolic Acid | N 2511AN | USP 2-0 45 cm 3/8 Circle Reverse Cutting 45 mm | 12 Foils |
| Polyglycolic Acid | N 2517 | USP 2-0 90 cm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2517N | USP 2-0 90 cm 3/8 Circle Reverse Cutting 45 mm | 12 Foils |
| Polyglycolic Acid - (U) | N 945 | USP 2-0 90 cm 1/2 Circle Round Body 36 mm | 12 Foils |
| Polyglycolic Acid | N 2382A | USP 2-0 45 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polyglycolic Acid | N 2473LDL | USP 2-0 100 cm 5/8 Circle Round Body 27 mm | 12 Foils |
| Polyglycolic Acid | N 2401 | USP 3-0 45 cm 3/8 Circle Cutting 22 mm | 12 Foils |
| Polyglycolic Acid - (U) | N 2407 | USP 3-0 70 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Polyglycolic Acid | N 2437 | USP 3-0 70 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Polyglycolic Acid | N 2471 | USP 3-0 45 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Polyglycolic Acid | N 2472 | USP 3-0 90 cm 1/2 Circle Cutting 22 mm | 12 Foils |
| Polyglycolic Acid | N 2495 | USP 3-0 45 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Polyglycolic Acid | N 2515 | USP 3-0 45 cm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Polyglycolic Acid | N 2516 | USP 3-0 90 cm 1/2 Circle Taper Cutting 36 mm | 12 Foils |
| Polyglycolic Acid | N 2437DA | USP 3-0 70 cm 3/8 Circle Round Body 22 mm 3/8 Circle Round Body 22 mm | 12 Foils |
| Polyglycolic Acid | N 2383 | USP 3-0 90 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglycolic Acid | N 2472VS | USP 3-0 110 cm 1/2 Circle Cutting 22 mm | 12 Foils |
| Polyglycolic Acid | N 2304 | USP 4-0 70 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Polyglycolic Acid | N 2465 | USP 4-0 45 cm 5/8 Circle Cutting 16 mm | 12 Foils |
| Polyglycolic Acid | N 2305 | USP 4-0 70 cm 1/2 Circle Round Body 16 mm | 12 Foils |
| Polyglycolic Acid | N 2494 | USP 4-0 45 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Polyglycolic Acid | N 2494N | USP 4-0 45 cm 3/8 Circle Reverse Cutting 16 mm | 12 Foils |
| Polyglycolic Acid | N 2944 | USP 4-0 70 cm 3/8 Circle Reverse Cutting 16 mm | 12 Foils |
| Polyglycolic Acid - (U) | N 2306 | USP 4-0 70 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Polyglycolic Acid | N 2304VS | USP 4-0 20 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Polyglycolic Acid | N 2303 | USP 5-0 45 cm 1/2 Circle Round Body 16 mm | 12 Foils |

| Product | Product Code | Description | Units / Box |
|-------------------|-----------------|---|-------------|
| Polyglycolic Acid | N 2442 | USP 5-0 45 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Polyglycolic Acid | N 2493 | USP 5-0 45 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Polyglycolic Acid | S 2615 | USP 1 150 cm Without Needle | 12 Foils |
| Polyglycolic Acid | S 2615MDT | USP 1 180 cm Without Needle | 12 Foils |
| Polyglycolic Acid | S 2614 | USP 0 150 cm Without Needle | 12 Foils |
| Polyglycolic Acid | S 2614MDT | USP 0 180 cm Without Needle | 12 Foils |
| Polyglycolic Acid | S 2613 | USP 2-0 150 cm Without Needle | 12 Foils |
| Polyglycolic Acid | S 2613MDT | USP 2-0 180 cm Without Needle | 12 Foils |
| Polyglycolic Acid | S 2777MDT/180 F | USP 2-0 180 cm Without Needle | 12 Foils |
| Polyglycolic Acid | S 2777MDT F | USP 2-0 150 cm Without Needle | 12 Foils |
| Polyglycolic Acid | N 2763 F | USP 0 90 cm 1/2 Circle Taper Cutting 35 mm | 12 Foils |
| Polyglycolic Acid | N 2763AECO | USP 0 70 cm 1/2 Circle Taper Cutting 35 mm | 12 Foils |
| Polyglycolic Acid | N 2763 NL F | USP 0 110 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglycolic Acid | N 2762 F | USP 2-0 90 cm 1/2 Circle Taper Cutting 35 mm | 12 Foils |
| Polyglycolic Acid | N 2762 FVS | USP 2-0 110 cm 1/2 Circle Taper Cutting 35 mm | 12 Foils |
| Polyglycolic Acid | N 2762 AF | USP 2-0 70 cm 1/2 Circle Taper Cutting 35 mm | 12 Foils |
| Polyglycolic Acid | N 2777 F | USP 2-0 140 cm 1/2 Circle Round Body 36 mm 1/2 Circle Reverse Cutting 36 mm | 12 Foils |
| Polyglycolic Acid | N 2777 F N | USP 2-0 140 cm 1/2 Circle Reverse Cutting 36 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyglycolic Acid | N 2762AS | USP 2-0 76 cm 1/2 Circle Reverse Cutting 36 mm | 12 Foils |
| Polyglycolic Acid | N 2732A F | USP 3-0 45 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Polyglycolic Acid | N 2732 F | USP 3-0 75 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Polyglycolic Acid | N 2732 F70 | USP 3-0 70 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Polyglycolic Acid | N 2735A F | USP 3-0 45 cm 1/2 Circle Cutting 22 mm | 12 Foils |

SanaGut

Catgut Chromic
Absorbable Surgical Suture U.S.P.
Sterile, Natural Surgical Suture



SanaGut is an absorbable, sterile surgical suture made from purified connective tissue (primarily collagen) derived from the submucosal fibrous layer of healthy sheep intestines. Treated with chromic salts, it offers enhanced resistance to enzymatic absorption, providing prolonged wound support. It ensures strong knot security, delayed breakdown, and superior tensile strength that lasts for 14–21 days, with complete natural absorption occurring in 90–120 days.

Applications include uterus closure, episiotomy repair, soft tissue approximation, and fast-healing tissues or muscles. It is a twisted, absorbable suture.

Specifications:

| | |
|--------------------|--|
| Type | : Twisted, Absorbable |
| Composition | : Catgut Chromic |
| Color | : Brown |
| Strength Retention | : 14–21 days |
| Absorption | : 90–120 days via enzymatic reaction |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 5 years |
| Packaging | : Aluminum foil packs (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

SanaGut Plain

Catgut Plain
Absorbable Surgical Suture U.S.P
Sterile, Natural Surgical Suture



SanaGut Plain is a fast-absorbing, sterile surgical suture made from a strand of collagenous material derived from the submucosal layers of healthy sheep intestines. It is designed for dermal suturing and external knot-tying procedures, eliciting minimal tissue reaction during absorption. It provides strong knot security, ideal tensile strength for 7–10 days, and complete natural absorption in 60–70 days.

Applications include uterus closure, episiotomy repair, soft tissue approximation, and fast-healing tissues or muscles. It is a twisted, absorbable suture.

Specifications:

| | |
|--------------------|--|
| Type | : Twisted, Absorbable |
| Composition | : Catgut Plain |
| Color | : Beige |
| Strength Retention | : 7–10 days |
| Absorption | : 60–70 days via enzymatic reaction |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 5 years |
| Packaging | : Aluminum foil packs (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

| Product | Product Code | Description | Units / Box |
|----------------|--------------|--|-------------|
| Catgut Chromic | N 4228 | USP 2 100 cm 1/2 Circle Round Body MF 45 mm | 12 Foils |
| Catgut Chromic | N 4244 | USP 2 100 cm 1/2 Circle Trocar Point (H) 40 mm | 12 Foils |
| Catgut Chromic | N 4278 | USP 2 100 cm 1/2 Circle Round Body 65 mm | 12 Foils |
| Catgut Chromic | N 4228A | USP 2 76 cm 1/2 Circle Round Body MF 45 mm | 12 Foils |
| Catgut Chromic | N 4228RC | USP 2 100 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Catgut Chromic | N 4226 | USP 1 76 cm 3/8 Circle Round Body 45 mm | 12 Foils |
| Catgut Chromic | N 4227 | USP 1 100 cm 1/2 Circle Round Body MF 45 mm | 12 Foils |
| Catgut Chromic | N 4239 | USP 1 100 cm 1/2 Circle Trocar Point (H) 40 mm | 12 Foils |
| Catgut Chromic | N 4239A | USP 1 76 cm 1/2 Circle Trocar Point (H) 40 mm | 12 Foils |
| Catgut Chromic | N 4259 | USP 1 76 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Catgut Chromic | N 4262 | USP 1 76 cm 3/8 Circle Cutting 60 mm | 12 Foils |
| Catgut Chromic | N 4267 | USP 1 76 cm 3/8 Circle Round Body Blunt 63 mm | 12 Foils |
| Catgut Chromic | N 4273 | USP 1 76 cm 1/2 Circle Reverse Cutting 50 mm | 12 Foils |
| Catgut Chromic | N 4275 | USP 1 152 cm 1/2 Circle Round Body 45 mm 1/2 Circle Round Body 45 mm | 12 Foils |
| Catgut Chromic | N 4297 | USP 1 100 cm 1/2 Circle Taper Cutting (H) 40 mm | 12 Foils |
| Catgut Chromic | N 4297A | USP 1 76 cm 1/2 Circle Taper Cutting (H) 40 mm | 12 Foils |
| Catgut Chromic | N 4826 | USP 1 76 cm 3/8 Circle Round Body 45 mm 3/8 Circle Round Body 45 mm | 12 Foils |
| Catgut Chromic | N 4226VS | USP 1 110 cm 3/8 Circle Round Body 45 mm | 12 Foils |
| Catgut Chromic | N 4114 | USP 0 76 cm 3/8 Circle Round Body 32 mm | 12 Foils |
| Catgut Chromic | N 4124 | USP 0 76 cm 1/2 Circle Round Body 26 mm | 12 Foils |
| Catgut Chromic | N 4205 | USP 0 76 cm Straight Round Body 45 mm | 12 Foils |
| Catgut Chromic | N 4213 | USP 0 76 cm Straight Round Body 60 mm | 12 Foils |
| Catgut Chromic | N 4217 | USP 0 76 cm 3/8 Circle Round Body 30 mm | 12 Foils |
| Catgut Chromic | N 4221 | USP 0 76 cm 3/8 Circle Round Body 40 mm | 12 Foils |
| Catgut Chromic | N 4225 | USP 0 76 cm 3/8 Circle Round Body 45 mm | 12 Foils |
| Catgut Chromic | N 4229 | USP 0 100 cm 1/2 Circle Round Body MF 45 mm | 12 Foils |
| Catgut Chromic | N 4233 | USP 0 76 cm 3/8 Circle Round Body 65 mm | 12 Foils |
| Catgut Chromic | N 4242 | USP 0 76 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Catgut Chromic | N 4246 | USP 0 76 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Catgut Chromic | N 4250 | USP 0 76 cm 1/2 Circle Round Body 45 mm | 12 Foils |

| Product | Product Code | Description | Units / Box |
|----------------|--------------|---|-------------|
| Catgut Chromic | N 4258 | USP 0 76 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Catgut Chromic | N 4264 | USP 0 76 cm 5/8 Circle Round Body 35 mm | 12 Foils |
| Catgut Chromic | N 4276 | USP 0 152 cm 1/2 Circle Round Body 45 mm 1/2 Circle Reverse Cutting 45 mm | 12 Foils |
| Catgut Chromic | N 4276N | USP 0 152 cm 1/2 Circle Round Body 40 mm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Catgut Chromic | N 4279 | USP 0 100 cm 1/2 Circle Taper Cutting (H) 40 mm | 12 Foils |
| Catgut Chromic | N 4296 | USP 0 76 cm 1/2 Circle Round Body 38 mm | 12 Foils |
| Catgut Chromic | N 4254 | USP 0 76 cm 1/2 Circle Round Body (H) 30 mm | 12 Foils |
| Catgut Chromic | N 4123 | USP 2-0 76 cm 1/2 Circle Round Body 26 mm | 12 Foils |
| Catgut Chromic | N 4200 | USP 2-0 76 cm Straight Round Body 40 mm | 12 Foils |
| Catgut Chromic | N 4204 | USP 2-0 76 cm Straight Round Body 45 mm | 12 Foils |
| Catgut Chromic | N 4208 | USP 2-0 76 cm Straight Round Body 50 mm | 12 Foils |
| Catgut Chromic | N 4216 | USP 2-0 76 cm 3/8 Circle Round Body 30 mm | 12 Foils |
| Catgut Chromic | N 4220 | USP 2-0 76 cm 3/8 Circle Round Body 40 mm | 12 Foils |
| Catgut Chromic | N 4232 | USP 2-0 76 cm 3/8 Circle Round Body 65 mm | 12 Foils |
| Catgut Chromic | N 4241 | USP 2-0 76 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Catgut Chromic | N 4245 | USP 2-0 76 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Catgut Chromic | N 4249 | USP 2-0 76 cm 1/2 Circle Round Body 45 mm | 12 Foils |
| Catgut Chromic | N 4255 | USP 2-0 76 cm 1/2 Circle Cutting (H) 30 mm | 12 Foils |
| Catgut Chromic | N 4277 | USP 2-0 152 cm 1/2 Circle Round Body 36 mm 1/2 Circle Reverse Cutting 36 mm | 12 Foils |
| Catgut Chromic | N 4295 | USP 2-0 76 cm 1/2 Circle Round Body 38 mm | 12 Foils |
| Catgut Plain | N 5033 | USP 2-0 76 cm 3/8 Circle Round Body 30 mm | 12 Foils |
| Catgut Plain | N 2003LDL | USP 2-0 152 cm 1/2 Circle Cutting 30 mm 1/2 Circle Cutting 30 mm | 12 Foils |
| Catgut Chromic | N 4201 | USP 3-0 76 cm 3/8 Circle Cutting 22 mm | 12 Foils |
| Catgut Chromic | N 4237 | USP 3-0 76 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Catgut Chromic | N 4272 | USP 3-0 76 cm 1/2 Circle Reverse Cutting 12 mm | 12 Foils |
| Catgut Chromic | N 4274 | USP 3-0 76 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Catgut Plain | N 4238 | USP 3-0 76 cm 1/2 Circle Round Body 25 mm | 12 Foils |
| Catgut Plain | N 4271 | USP 3-0 76 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Catgut Plain | N 4288 | USP 3-0 76 cm Straight Cutting 60 mm | 12 Foils |
| Catgut Plain | N 5007 | USP 3-0 76 cm Straight Round Body 50 mm | 12 Foils |

| Product | Product Code | Description | Units / Box |
|----------------|--------------|--|-------------|
| Catgut Plain | N 5010 | USP 3-0 76 cm 3/8 Circle Cutting 22 mm | 12 Foils |
| Catgut Plain | N 5020 | USP 3-0 76 cm 3/8 Circle Cutting 26 mm | 12 Foils |
| Catgut Chromic | N 4223 | USP 3-0 76 cm 3/8 Circle Round Body 45 mm | 12 Foils |
| Catgut Chromic | N 4265 | USP 3-0 38 cm 5/8 Circle Cutting 16 mm | 12 Foils |
| Catgut Chromic | N 4280 | USP 4-0 76 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Catgut Chromic | N 5048 | USP 4-0 76 cm 3/8 Circle Round Body 16 mm | 12 Foils |
| Catgut Plain | N 4270 | USP 4-0 76 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Catgut Chromic | N 4266 | USP 4-0 38 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Catgut Chromic | N 4268 | USP 5-0 76 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Catgut Chromic | N 4281 | USP 5-0 76 cm 3/8 Circle Round Body 16 mm | 12 Foils |
| Catgut Plain | N 4269 | USP 5-0 76 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Catgut Plain | S 2006 | USP 2 152 cm Without Needle | 12 Foils |
| Catgut Plain | S 2005 | USP 1 152 cm Without Needle | 12 Foils |
| Catgut Plain | S 2004 | USP 0 152 cm Without Needle | 12 Foils |
| Catgut Plain | S 2003 | USP 2-0 152 cm Without Needle | 12 Foils |
| Catgut Plain | S 2002 | USP 3-0 152 cm Without Needle | 12 Foils |
| Catgut Chromic | S 2242 | USP 3-0 150 cm Without Needle | 12 Foils |
| Catgut Chromic | S 2217 | USP 3 152 cm Without Needle | 12 Foils |
| Catgut Chromic | S 2216 | USP 2 152 cm Without Needle | 12 Foils |
| Catgut Chromic | S 2215 | USP 1 152 cm Without Needle | 12 Foils |
| Catgut Chromic | S 2255 | USP 1 150 cm Without Needle | 12 Foils |
| Catgut Chromic | S 2214 | USP 0 152 cm Without Needle | 12 Foils |
| Catgut Chromic | S 2213 | USP 2-0 152 cm Without Needle | 12 Foils |
| Catgut Chromic | S 2212 | USP 3-0 152 cm Without Needle | 12 Foils |

SanaCryl Mono

Monofilament Polyglactaprone
Absorbable Surgical Suture
Sterile, Synthetic Surgical Suture



SanaCryl Mono is a synthetic, absorbable, sterile monofilament suture made from a copolymer of 75% glycolide and 25% epsilon-caprolactone. It is designed for general soft tissue approximation and ligation, making it ideal for female skin closure and subcutaneous approximation. However, it is not recommended for cardiovascular, neurological, Sanasurgical, or ophthalmic procedures. It offers strong knot security, a smooth surface for minimal tissue reaction, and superior tensile strength that retains 50% strength for 14 days, with full absorption in 90–120 days.

Applications include female skin approximation and subcutaneous approximation. It is a monofilament, synthetic absorbable suture.

Specifications:

| | |
|--------------------|--|
| Type | : Monofilament, Synthetic Absorbable |
| Composition | : 75% Poly(glycolide) – 25% Co-Caprolactone |
| Coating | : None |
| Color | : Undyed / Violet |
| Strength Retention | : 50% (14 days), provides wound support for 21–28 days |
| Absorption | : 90–120 days via hydrolysis |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 3 years |
| Packaging | : Aluminum foil packs (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

| Product | Product Code | Description | Units / Box |
|-----------------------|--------------|---|-------------|
| Polyglecaprone 25 (U) | S 1614 | USP 0 150 cm Without Needle | 12 Foils |
| Polyglecaprone 25 (U) | N 1642 | USP 0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglecaprone 25 (U) | N 1641 | USP 2-0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglecaprone 25 (U) | N 1651/35 | USP 2-0 35 cm Straight Cutting 60 mm | 12 Foils |
| Polyglecaprone 25 (U) | N 1741 | USP 2-0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polyglecaprone 25 (U) | N 1326 | USP 3-0 70 cm 3/8 Circle Precision Cutting MF 25 mm | 12 Foils |
| Polyglecaprone 25 (U) | N 1326A | USP 3-0 45 cm 3/8 Circle Precision Cutting MF 25 mm | 12 Foils |
| Polyglecaprone 25 (U) | N 3326PS | USP 3-0 70 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Polyglecaprone 25 (U) | N 1650 | USP 3-0 70 cm Straight Cutting 60 mm | 12 Foils |
| Polyglecaprone 25 (U) | N 1201 | USP 3-0 90 cm 1/2 Circle Round Body 35 mm | 12 Foils |
| Polyglecaprone 25 (U) | N 1205 | USP 4-0 70 cm 3/8 Circle Reverse Cutting 16 mm | 12 Foils |

SanaSynth PDO

Monofilament Polydioxanone
Absorbable Surgical Suture
Sterile, Synthetic Surgical Suture



SanaSynth PDO is a synthetic, absorbable, sterile monofilament suture made from polyester poly(p-dioxanone). Designed for extended wound support, it is ideal for procedures requiring delayed healing, including surgeries on geriatric, obese, diabetic, and oncology patients. It provides delayed wound support, minimal tissue drag due to its smooth surface, and superior tensile strength that retains 50% of its strength for 30 days, with full absorption occurring in 180–210 days.

Applications include joint muscle approximation, hard tissue approximation, and gluteus maximus approximation. It is a monofilament, synthetic absorbable suture.

Specifications:

| | |
|--------------------|--|
| Type | : Monofilament, Synthetic Absorbable |
| Composition | : Polyester, Poly(p-dioxanone) |
| Coating | : None |
| Color | : Violet |
| Strength Retention | : 50% (30 days), provides wound support for 35–45 days |
| Absorption | : 180–210 days via hydrolysis |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 3 years |
| Packaging | : Aluminum foil packs (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

| Product | Product Code | Description | Units / Box |
|----------------------|--------------|---|-------------|
| Polydioxanone | N 9234 | USP 1 90 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polydioxanone | N 9248 | USP 1 90 cm 1/2 Circle Round Body 45 mm | 12 Foils |
| Polydioxanone | N 9352 | USP 1 90 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polydioxanone | N 9352PS | USP 1 45 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polydioxanone | N 9221 | USP 1 90 cm 1/2 Circle Reverse Cutting (H) 40 mm | 12 Foils |
| Polydioxanone | N 9235SA | USP 1 110 cm 1/2 Circle Taper Cutting 30 mm | 12 Foils |
| Polydioxanone | N 9826 | USP 1 70 cm 1/2 Circle Reverse Cutting 23 mm | 12 Foils |
| Polydioxanone | N 9235 | USP 1 90 cm 1/2 Circle Taper Cutting 40 mm | 12 Foils |
| Polydioxanone | N 9262 | USP 1 LOOP 150 cm 1/2 Circle Round Body (H) 50 mm | 12 Foils |
| Polydioxanone | N 9210 | USP 0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polydioxanone | N 9233 | USP 0 70 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polydioxanone | N 9371 | USP 0 90 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polydioxanone | N 9261 | USP 0 LOOP 150 cm 1/2 Circle Round Body (H) 50 mm | 12 Foils |
| Polydioxanone | N 9133 | USP 2-0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polydioxanone | N 9132 | USP 3-0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polydioxanone | N 9237 | USP 3-0 70 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Polydioxanone | N 9238 | USP 3-0 70 cm 3/8 Circle Reverse Cutting 25 mm | 12 Foils |
| Polydioxanone | N 9304 | USP 4-0 70 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Polydioxanone | N 9109 | USP 4-0 90 cm 1/2 Circle Round Body 17 mm 1/2 Circle Round Body 17 mm | 12 Foils |
| Polydioxanone | N 9201H | USP 5-0 70 cm 1/2 Circle Round Body 13 mm 1/2 Circle Round Body 13 mm | 12 Foils |
| Polydioxanone | N 9201 | USP 5-0 70 cm 3/8 Circle Round Body 13 mm | 12 Foils |
| Polydioxanone | N 9091 | USP 6-0 75 cm 3/8 Circle Round Body 8 mm 3/8 Circle Round Body 8 mm | 12 Foils |
| Barbed Polydioxanone | BS N0315U | USP 2-0 30 cm 1/2 Circle Round Body 37 mm | 6 Foils |
| Barbed Polydioxanone | BS N0316U | USP 0 30 cm 1/2 Circle Round Body 37 mm | 6 Foils |
| Barbed Polydioxanone | BS N0325U | USP 2-0 45 cm 1/2 Circle Round Body 37 mm | 6 Foils |
| Barbed Polydioxanone | BS N0326U | USP 0 45 cm 1/2 Circle Round Body 37 mm | 6 Foils |
| Barbed Polydioxanone | BS N0615U | USP 2-0 30 cm 1/2 Circle Round Body 26 mm | 6 Foils |
| Barbed Polydioxanone | BS N0624U | USP 3-0 45 cm 1/2 Circle Round Body 26 mm | 6 Foils |
| Barbed Polydioxanone | BS N1100U | USP 3-0 45 cm 3/8 Circle Reverse Cutting 26 mm | 6 Foils |

SanaLene

Monofilament Polypropylene
Non-Absorbable Surgical Suture U.S.P.
Sterile, Synthetic Surgical Suture



SanaLene is a non-absorbable, sterile surgical suture made from isotactic crystalline stereo isomer of polypropylene, a synthetic linear polyolefin. Pigmented blue for enhanced visibility, it is designed for general soft tissue approximation and ligation, including cardiovascular, ophthalmic, and neurological procedures. It features a smooth surface for minimal tissue trauma, monofilament structure to prevent bacterial growth and capillary action, uniform elasticity that supports movement for better cosmetic outcomes, and enhanced visibility due to its blue pigmentation.

Applications include skin closure, hard tissue approximation, ligament repair, and mesh fixation. It is a monofilament, non-absorbable suture.

Specifications:

| | |
|---------------|--|
| Type | : Monofilament, Non-Absorbable |
| Composition | : Polypropylene |
| Color | : Blue |
| Absorption | : Non-Absorbable |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 5 years |
| Packaging | : Food-grade folder packed in a medical-grade pouch (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

| Product | Product Code | Description | Units / Box |
|---------------|--------------|---|-------------|
| Polypropylene | N 840 | USP 1 70 cm 1/2 Circle Round Body (H) 30 mm | 12 Foils |
| Polypropylene | N 843 | USP 1 70 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polypropylene | N 843A | USP 1 35 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polypropylene | N 843VS | USP 1 90 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polypropylene | N 845 | USP 1 100 cm 1/2 Circle Taper Cutting (H) 45 mm | 12 Foils |
| Polypropylene | N 883 | USP 1 100 cm 1/2 Circle Reverse Cutting (H) 45 mm | 12 Foils |
| Polypropylene | N 896 | USP 1 100 cm 1/2 Circle Taper Cutting (H) 30 mm | 12 Foils |
| Polypropylene | N 883MB | USP 1 45 cm 1/2 Circle Reverse Cutting (H) 40 mm | 12 Foils |
| Polypropylene | N 883DB | USP 1 45 cm 1/2 Circle Reverse Cutting 40 mm | 12 Foils |
| Polypropylene | N 805 | USP 0 70 cm 1/2 Circle Taper Cutting 25 mm | 12 Foils |
| Polypropylene | N 830 | USP 0 70 cm 1/2 Circle Round Body (H) 30 mm | 12 Foils |
| Polypropylene | N 842 | USP 0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polypropylene | N 842A | USP 0 35 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polypropylene | N 842VS | USP 0 90 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polypropylene | N 894 | USP 0 90 cm 1/2 Circle Taper Cutting (H) 30 mm | 12 Foils |
| Polypropylene | N 824 | USP 0 90 cm 1/2 Circle Round Body 26 mm 1/2 Circle Round Body 26 mm | 12 Foils |
| Polypropylene | N 806 | USP 2-0 90 cm 1/2 Circle Taper Cutting 25 mm | 12 Foils |
| Polypropylene | N 807 | USP 2-0 70 cm 1/2 Circle Taper Cutting 25 mm | 12 Foils |
| Polypropylene | N 841 | USP 2-0 70 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polypropylene | N 841VS | USP 2-0 90 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Polypropylene | N 898 | USP 2-0 70 cm 1/2 Circle Taper Cutting 17 mm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Polypropylene | N 018 | USP 3-0 70 cm 3/8 Circle Cutting 22 mm | 12 Foils |
| Polypropylene | N 884 | USP 2-0 70 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Polypropylene | N 822 | USP 3-0 70 cm 1/2 Circle Round Body 25 mm 1/2 Circle Round Body 25 mm | 12 Foils |
| Polypropylene | N 825 | USP 3-0 90 cm 1/2 Circle Round Body 25 mm | 12 Foils |
| Polypropylene | N 838 | USP 3-0 70 cm 1/2 Circle Round Body 25 mm | 12 Foils |
| Polypropylene | N 887 | USP 3-0 70 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Polypropylene | N 888 | USP 3-0 70 cm Straight Cutting 60 mm | 12 Foils |
| Polypropylene | N 849 | USP 4-0 70 cm 1/2 Circle Round Body 16 mm | 12 Foils |
| Polypropylene | N 853 | USP 4-0 90 cm 1/2 Circle Taper Cutting 17 mm 1/2 Circle Taper Cutting 17 mm | 12 Foils |

| Product | Product Code | Description | Units / Box |
|---------------|--------------|---|-------------|
| Polypropylene | N 867 | USP 4-0 90 cm 1/2 Circle Taper Cutting 25 mm 1/2 Circle Taper Cutting 25 mm | 12 Foils |
| Polypropylene | N 870PS | USP 4-0 70 cm 3/8 Circle Cutting MicroPrime 16 mm | 12 Foils |
| Polypropylene | N 879 | USP 4-0 45 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Polypropylene | N 857 | USP 4-0 90 cm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Polypropylene | N 880 | USP 5-0 70 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Polypropylene | N 881 | USP 5-0 70 cm 3/8 Circle Round Body 16 mm | 12 Foils |
| Polypropylene | N 882 | USP 5-0 70 cm 3/8 Circle Round Body 16 mm 3/8 Circle Round Body 16 mm | 12 Foils |
| Polypropylene | N 880VS | USP 5-0 90 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Polypropylene | N 826 | USP 5-0 70 cm 3/8 Circle Slim Blade Cutting 15 mm | 12 Foils |
| Polypropylene | N 892 | USP 5-0 90 cm 3/8 Circle Round Body 20 mm 3/8 Circle Round Body 20 mm | 12 Foils |
| Polypropylene | N 823 | USP 6-0 70 cm 3/8 Circle Slim Blade Cutting 15 mm | 12 Foils |
| Polypropylene | N 829 | USP 6-0 70 cm 3/8 Circle Round Body 13 mm 3/8 Circle Round Body 13 mm | 12 Foils |
| Polypropylene | N 878 | USP 6-0 45 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Polypropylene | N 878VS | USP 6-0 70 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |

SanaBond

Braided Polyester
Non-Absorbable Surgical Suture U.S.P.
Sterile, Synthetic Surgical Suture



SanaBond is a non-absorbable, braided, and sterile surgical suture made from high-molecular-weight poly(ethylene terephthalate). Designed for optimal handling and high visibility in the surgical field, it provides superior strength and durability for a range of soft tissue approximation and ligation procedures, including cardiovascular, ophthalmic, and neurological surgeries. It features an ultra-strong fiber for enhanced support, a smooth surface for easy passage through tough tissues, and excellent handling characteristics preferred by orthopedic surgeons.

Applications include general soft tissue approximation, cardiovascular procedures, and orthopedic surgeries. It is a braided, non-absorbable suture.

Specifications:

| | |
|---------------|--|
| Type | : Braided, Non-Absorbable |
| Composition | : Polyester |
| Color | : Green / White |
| Absorption | : Non-Absorbable |
| USP Range | : 5/0 – 5 |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 5 years |
| Packaging | : Food-grade folder packed in a medical-grade pouch (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

| Product | Product Code | Description | Units / Box |
|-------------------|--------------|---|-------------|
| Polyester (Green) | N 645 | USP 5 100 cm 1/2 Circle Trocar Point (H) 55 mm | 12 Foils |
| Polyester (Green) | N 643 | USP 2 100 cm 1/2 Circle Taper Cutting (H) 45 mm | 12 Foils |
| Polyester (Green) | N 680 | USP 1 100 cm 1/2 Circle Reverse Cutting 50 mm | 12 Foils |
| Polyester (Green) | N 679 | USP 0 100 cm 1/2 Circle Reverse Cutting 50 mm | 12 Foils |
| Polyester (Green) | N 617 | USP 2-0 90 cm 1/2 Circle Taper Cutting 17 mm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Polyester (White) | N 617A | USP 2-0 90 cm 1/2 Circle Taper Cutting 17 mm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Polyester (Green) | N 637 | USP 2-0 76 cm 1/2 Circle Taper Cutting 17 mm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Polyester (Green) | N 677 | USP 2-0 76 cm 1/2 Circle Taper Cutting 25 mm 1/2 Circle Taper Cutting 25 mm | 12 Foils |
| Polyester (White) | N 677A | USP 2-0 76 cm 1/2 Circle Taper Cutting 25 mm 1/2 Circle Taper Cutting 25 mm | 12 Foils |
| Polyester (Green) | N 687 | USP 2-0 90 cm 1/2 Circle Taper Cutting 25 mm 1/2 Circle Taper Cutting 25 mm | 12 Foils |
| Polyester (White) | N 687A | USP 2-0 90 cm 1/2 Circle Taper Cutting 25 mm 1/2 Circle Taper Cutting 25 mm | 12 Foils |
| Polyester (White) | N 687 | USP 2-0 75 cm 1/2 Circle Taper Cutting 26 mm 1/2 Circle Taper Cutting 26 mm | 12 Foils |
| Polyester (Green) | N 625 | USP 3-0 90 cm 1/2 Circle Round Body 25 mm 1/2 Circle Round Body 25 mm | 12 Foils |
| Polyester (Green) | N 632 | USP 3-0 76 cm 1/2 Circle Round Body 25 mm | 12 Foils |
| Polyester (Green) | N 636 | USP 3-0 76 cm 1/2 Circle Taper Cutting 17 mm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Polyester (Green) | N 652 | USP 3-0 76 cm 1/2 Circle Round Body 25 mm 1/2 Circle Round Body 25 mm | 12 Foils |
| Polyester (Green) | N 669 | USP 3-0 90 cm 1/2 Circle Taper Cutting 17 mm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Polyester (Green) | N 635 | USP 4-0 76 cm 1/2 Circle Taper Cutting 17 mm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Polyester (Green) | N 640 | USP 4-0 45 cm Straight Cutting 53 mm Straight Cutting 53 mm | 12 Foils |
| Polyester (Green) | N 672 | USP 4-0 76 cm 3/8 Circle Round Body 16 mm | 12 Foils |
| Polyester (Green) | N 691 | USP 4-0 90 cm 3/8 Circle Round Body 12 mm 3/8 Circle Round Body 12 mm | 12 Foils |
| Polyester (Green) | N 695 | USP 4-0 90 cm 1/2 Circle Taper Cutting 17 mm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Polyester (Green) | N 690 | USP 5-0 76 cm 3/8 Circle Round Body 12 mm 3/8 Circle Round Body 12 mm | 12 Foils |
| Polyester (Green) | N 696 | USP 5-0 76 cm 3/8 Circle Round Body 20 mm 3/8 Circle Round Body 20 mm | 12 Foils |

SanaLon

Monofilament Polyamide
Non-Absorbable Surgical Suture U.S.P.
Sterile, Synthetic Surgical Suture



SanaLon is a non-absorbable, sterile surgical suture made from long-chain aliphatic homopolymers of Nylon 6 and Nylon 6,6. Dyed black or blue for enhanced visibility, it is ideal for general soft tissue approximation and ligation, including cardiovascular, ophthalmic, and neurological procedures. It offers strong knot security, a smooth surface that prevents bacterial growth and capillary action, uniform elasticity for better cosmetic outcomes, and an ultra-smooth finish to reduce patient discomfort during suturing.

Applications include skin closure, hard tissue approximation, and casualty stitching. It is a monofilament, non-absorbable suture.

Specifications:

| | |
|---------------|--|
| Type | : Monofilament, Non-Absorbable |
| Composition | : Polyamide (Nylon) |
| Color | : Black or Blue |
| Absorption | : Non-Absorbable |
| USP Range | : 6/0 – 2 |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 5 years |
| Packaging | : Food-grade folder packed in a medical-grade pouch (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

| Product | Product Code | Description | Units / Box |
|-----------|--------------|---|-------------|
| Polyamide | N 3306 | USP 2 75 cm 1/2 Circle Reverse Cutting 38 mm | 12 Foils |
| Polyamide | N 3338 | USP 1 100 cm 1/2 Circle Reverse Cutting (H) 45 mm | 12 Foils |
| Polyamide | N 3338A | USP 1 45 cm 1/2 Circle Reverse Cutting (H) 45 mm | 12 Foils |
| Polyamide | N 3347 | USP 1 100 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyamide | N 3305 | USP 1 75 cm 1/2 Circle Reverse Cutting 38 mm | 12 Foils |
| Polyamide | N 3348 Loop | USP 1 150 cm 1/2 Circle Round Body (H) 50 mm | 12 Foils |
| Polyamide | N 3304 | USP 0 75 cm 1/2 Circle Reverse Cutting 38 mm | 12 Foils |
| Polyamide | N 3346 | USP 0 100 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Polyamide | N 3340 Loop | USP 0 150 cm 1/2 Circle Round Body 40 mm | 12 Foils |
| Polyamide | N 3389 | USP 0 70 cm Straight Cutting 60 mm | 12 Foils |
| Polyamide | N 3337 | USP 0 70 cm 3/8 Circle Reverse Cutting 45 mm | 12 Foils |
| Polyamide | N 3303 | USP 2-0 75 cm 1/2 Circle Reverse Cutting 38 mm | 12 Foils |
| Polyamide | N 3336 | USP 2-0 70 cm 3/8 Circle Reverse Cutting 45 mm | 12 Foils |
| Polyamide | N 3336A | USP 2-0 35 cm 3/8 Circle Reverse Cutting 45 mm | 12 Foils |
| Polyamide | N 3336DK | USP 2-0 110 cm 3/8 Circle Reverse Cutting 45 mm | 12 Foils |
| Polyamide | N 3390 | USP 2-0 70 cm Straight Cutting 60 mm | 12 Foils |
| Polyamide | N 3321 | USP 3-0 70 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Polyamide | N 3321A | USP 3-0 35 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Polyamide | N 3328 | USP 3-0 70 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Polyamide | N 3328A | USP 3-0 35 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Polyamide | N 3328DK | USP 3-0 110 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Polyamide | N 3388 | USP 3-0 70 cm Straight Cutting 60 mm | 12 Foils |
| Polyamide | N 3318 | USP 4-0 70 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Polyamide | N 3326N | USP 4-0 38 cm 3/8 Circle Cutting 10 mm | 12 Foils |
| Polyamide | N 3318S | USP 4-0 35 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Polyamide | N 3318RC | USP 4-0 70 cm 3/8 Circle Reverse Cutting 16 mm | 12 Foils |
| Polyamide | N 3319PS | USP 4-0 90 cm 3/8 Circle Cutting MicroPrime 16 mm | 12 Foils |
| Polyamide | N 3315 | USP 5-0 35 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Polyamide | N 3317 | USP 5-0 70 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Polyamide | N 3331 | USP 6-0 70 cm 3/8 Circle Cutting 26 mm | 12 Foils |

| Product | Product Code | Description | Units / Box |
|-----------|--------------|---|-------------|
| Polyamide | N 3719 | USP 10-0 38 cm 3/8 Circle Micro Point Spatulated 6 mm X 2 Nos | 12 Foils |
| Polyamide | S 903 | USP 2-0 150 cm Without Needle | 12 Foils |
| Polyamide | S 904 | USP 0 150 cm Without Needle | 12 Foils |
| Polyamide | S 905 | USP 1 150 cm Without Needle | 12 Foils |
| Polyamide | S 906 | USP 2 150 cm Without Needle | 12 Foils |

SanaSil

Black Braided Silk
Non-Absorbable Surgical Suture U.S.P.
Sterile, Natural Surgical Suture



SanaSil is a non-absorbable, sterile surgical suture made from fibroin, an organic protein derived from the domesticated silkworm (*Bombyx mori*). Processed to remove natural waxes and gums, it is coated with a special wax mixture to ensure smooth handling. Available in black and natural colors, it delivers excellent knot security and superior handling properties. It features high-purity protein to minimize tissue reaction, a smooth wax coating that reduces suturing time and trauma, and enhanced visibility with USP-approved dye.

Applications include skin closure, hard tissue approximation, tendon repair, and muscle approximation. It is a braided, non-absorbable suture.

Specifications:

| | |
|---------------|--|
| Type | : Braided, Non-Absorbable |
| Composition | : Braided Silk |
| Color | : Black (also available in natural) |
| Absorption | : Non-Absorbable |
| USP Range | : 5/0 – 2 |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 5 years |
| Packaging | : Food-grade folder packed in a medical-grade pouch (with/without needles) |
| Needle Types | : Round-bodied, cutting edge, taper point, straight, blunt point |

| Product | Product Code | Description | Units / Box |
|--------------------|--------------|--|-------------|
| Black Braided Silk | N 5061Pcm | USP 4 76 cm 1/2 Circle Round Body 35 mm | 12 Foils |
| Black Braided Silk | N 5061NN | USP 4 45 cm 1/2 Circle Taper Cutting 35 mm | 12 Foils |
| Black Braided Silk | N 5006 | USP 2 75 cm 1/2 Circle Reverse Cutting 38 mm | 12 Foils |
| Black Braided Silk | N 5006A | USP 2 76 cm 1/2 Circle Reverse Cutting 65 mm | 12 Foils |
| Black Braided Silk | N 5005 | USP 1 75 cm 1/2 Circle Reverse Cutting 38 mm | 12 Foils |
| Black Braided Silk | N 5062 | USP 1 76 cm 3/8 Circle Cutting 60 mm | 12 Foils |
| Black Braided Silk | N 5063 | USP 1 76 cm 1/2 Circle Round Body (H) 40 mm | 12 Foils |
| Black Braided Silk | N 5004 | USP 0 75 cm 1/2 Circle Reverse Cutting 38 mm | 12 Foils |
| Black Braided Silk | N 5017 | USP 0 76 cm 3/8 Circle Round Body 30 mm | 12 Foils |
| Black Braided Silk | N 5037 | USP 0 76 cm 3/8 Circle Reverse Cutting 45 mm | 12 Foils |
| Black Braided Silk | N 5037A | USP 0 35 cm 3/8 Circle Reverse Cutting 45 mm | 12 Foils |
| Black Braided Silk | N 5334 | USP 0 76 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Black Braided Silk | N 5333VS | USP 2-0 90 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Black Braided Silk | N 5335 | USP 0 75 cm 3/8 Circle Round Body 25 mm | 12 Foils |
| Black Braided Silk | N 5003 | USP 2-0 75 cm 1/2 Circle Reverse Cutting 38 mm | 12 Foils |
| Black Braided Silk | N 5036 | USP 2-0 76 cm 3/8 Circle Reverse Cutting 45 mm | 12 Foils |
| Black Braided Silk | N 5036A | USP 2-0 35 cm 3/8 Circle Reverse Cutting 45 mm | 12 Foils |
| Black Braided Silk | N 5052 | USP 2-0 76 cm 1/2 Circle Reverse Cutting 50 mm | 12 Foils |
| Black Braided Silk | N 5290 | USP 2-0 76 cm Straight Cutting 60 mm | 12 Foils |
| Black Braided Silk | N 5332 | USP 2-0 76 cm 1/2 Circle Round Body 26 mm | 12 Foils |
| Black Braided Silk | N 5333 | USP 2-0 76 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Black Braided Silk | N 5333KG | USP 2-0 90 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Black Braided Silk | N 5610 | USP 2-0 76 cm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Black Braided Silk | N 5670 | USP 2-0 76 cm 1/2 Circle Taper Cutting 25 mm | 12 Foils |
| Black Braided Silk | N 5002 | USP 3-0 76 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Black Braided Silk | N 5002VS | USP 3-0 90 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Black Braided Silk | N 5028 | USP 3-0 76 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Black Braided Silk | N 5028A | USP 3-0 35 cm 3/8 Circle Reverse Cutting 26 mm | 12 Foils |
| Black Braided Silk | N 5070 | USP 3-0 76 cm 1/2 Circle Round Body 25 mm | 12 Foils |
| Black Braided Silk | N 5083 | USP 3-0 76 cm 1/2 Circle Round Body 35 mm | 12 Foils |

| Product | Product Code | Description | Units / Box |
|--------------------------|--------------|--|-------------|
| Black Braided Silk | N 5087 | USP 3-0 76 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Black Braided Silk | N 5085 | USP 3-0 90 cm 1/2 Circle Round Body 25 mm | 12 Foils |
| Black Braided Silk | N 5332N | USP 3-0 76 cm 1/2 Circle Round Body 30 mm | 12 Foils |
| Black Braided Silk | N 5001 | USP 4-0 76 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Black Braided Silk | N 5049 | USP 4-0 76 cm 1/2 Circle Round Body 16 mm | 12 Foils |
| Black Braided Silk | N 5082 | USP 4-0 76 cm 3/8 Circle Round Body 16 mm | 12 Foils |
| Black Braided Silk | N 5086 | USP 4-0 76 cm 1/2 Circle Round Body 20 mm | 12 Foils |
| Black Braided Silk | N 5001VS | USP 4-0 90 cm 3/8 Circle Cutting 16 mm | 12 Foils |
| Black Braided Silk | N 5080 | USP 5-0 76 cm 3/8 Circle Reverse Cutting 12 mm | 12 Foils |
| Black Braided Silk | N 5081 | USP 5-0 76 cm 3/8 Circle Round Body 16 mm | 12 Foils |
| Black Braided Silk | N 5090 | USP 5-0 76 cm Straight Round Body 20 mm | 12 Foils |
| Black Braided Silk | N 5096 | USP 5-0 76 cm 3/8 Circle Round Body 20 mm | 12 Foils |
| Black Braided Silk | N 5098 | USP 5-0 76 cm 1/2 Circle Taper Cutting 17 mm | 12 Foils |
| Black Braided Silk | N 5027 | USP 5-0 76 cm 3/8 Circle Round Body 12 mm | 12 Foils |
| Black Braided Silk | N 5029 | USP 6-0 76 cm 3/8 Circle Round Body 13 mm | 12 Foils |
| Black Braided Silk | S 211 | USP 4-0 150 cm Without Needle | 12 Foils |
| Black Braided Silk | S 212 | USP 3-0 150 cm Without Needle | 12 Foils |
| Black Braided Silk | S 213 | USP 2-0 150 cm Without Needle | 12 Foils |
| Black Braided Silk | S 214 | USP 0 150 cm Without Needle | 12 Foils |
| Black Braided Silk | S 215 | USP 1 150 cm Without Needle | 12 Foils |
| Black Braided Silk | S 216 | USP 2 150 cm Without Needle | 12 Foils |
| Black Braided Silk | S 217 | USP 3 150 cm Without Needle | 12 Foils |
| Black Braided Silk Reels | R 826 | USP 2 2500 cm | 6 Reels |
| Black Braided Silk Reels | R 825 | USP 1 2500 cm | 6 Reels |
| Black Braided Silk Reels | R 824 | USP 0 2500 cm | 6 Reels |
| Black Braided Silk Reels | R 823 | USP 2-0 2500 cm | 6 Reels |
| Black Braided Silk Reels | R 822 | USP 3-0 2500 cm | 6 Reels |
| Black Braided Silk Reels | R 821 | USP 4-0 2500 cm | 6 Reels |

SanaLene Mesh

Monofilament Polypropylene Mesh for
Hernia Repair Sterile, Synthetic,
Non-Absorbable Surgical Mesh



SanaLene Mesh is a sterile, non-absorbable surgical mesh made from knitted filaments of extruded polypropylene, identical in composition to SanaLene Polypropylene Suture. It is specifically designed for hernia repair and other fascial deficiencies that require reinforcement. Its unique bi-directional elasticity allows it to adapt naturally to body movements, and it can be cut into any shape or size without unraveling. It offers strong and durable support, is flexible and lightweight for enhanced patient comfort, and features an ultra-smooth surface to reduce tissue irritation and promote healing.

Applications include hernia repair and reinforcement of traumatic or surgical wounds. It is a non-absorbable polypropylene mesh.

Specifications:

| | |
|----------------|--|
| Type | : Knitted Filaments of Extruded Polypropylene |
| Composition | : Polypropylene |
| Color | : White (Undyed) |
| Absorption | : Non-Absorbable |
| Sterilization | : Ethylene Oxide (EO) |
| Shelf Life | : 5 years |
| Standard Sizes | : 6 x 11 cm, 8 x 15 cm, 10 x 15 cm, 15 x 15 cm, 30 x 30 cm |
| Customization | : Any size and shape available upon request |

SanaWax

Sterile BoneWax 2.5g



SanaWax is a sterile bone wax composed of beeswax, paraffin, and isopropyl palmitate, a wax-softening agent. It is used to control bleeding from bone surfaces during surgical procedures by forming a mechanical barrier (tamponade). SanaWax is non-absorbable, biochemically inert, and remains at the application site to aid in local haemostasis. It provides effective haemostasis by blocking bleeding, is inert so it doesn't interfere with healing, and comes sterile and ready to use in individual foil packs to prevent contamination.

Applications include control of bone surface bleeding during surgical, orthopedic, and neurosurgical procedures. It is a non-absorbable bone haemostatic agent.

Specifications:

| | |
|--------------------|---|
| Material | : Beeswax, Paraffin & Isopropyl Palmitate |
| Color | : Opaque |
| Absorption | : Minimally Resorbable |
| Sterilization | : Gamma Radiation |
| Shelf Life | : 5 years |
| Packaging | : Aluminum foil packs (2.5g per pack) |
| Storage Conditions | : Store in a cool, dry place away from excessive heat |

| Product | Product Code | Description | Units / Box |
|-----------------|--------------|--|---------------|
| Hernia Mesh | MM 0611 | 6 cm x 11 cm | 6 Foils / Box |
| Hernia Mesh | MM 0715 | 7.6 cm x 15 cm | 6 Foils / Box |
| Hernia Mesh | MM 1015 | 10 cm x 15 cm | 3 Foils / Box |
| Hernia Mesh | MM 1515 | 15 cm x 15 cm | 3 Foils / Box |
| Hernia Mesh | MM 3030 | 30 cm x 30 cm | 1 Foil / Box |
| Hernia Kit | MK 611 | MM 6cm X 11cm MM 0611, N 018, N 841 & N 1326 | 6 Kits / Box |
| Hernia Kit | MK 715 | MM 7.6cm X 15cm MM 0715, N 018, N 841 & N 1326 | 6 Kits / Box |
| Catgut Endoloop | ML 2215 | Pre-Knotted Loop Ligature | 6 Foils / Box |
| Bone Wax | S 810 | Sterile Bone Wax | 12 Foils |

Our Products

| | | |
|-------------------|------------------------------|--------------------------|
| 1. SanaSynth | - Braided and Coated | - Polyglactin 910 |
| 2. SanaSynth Fast | - Braided and Coated | - Polyglactin 910 Fast |
| 3. SanaCryl | - Braided and Coated | - Polyglycolic Acid |
| 4. SanaCryl Fast | - Braided and Coated | - Polyglycolic Acid Fast |
| 5. SanaSynth PDO | - Monofilamen | - Polydioxanone |
| 6. SanaCryl Mono | - Monofilament | - Poliglecaprone |
| 7. SanaGut | - Twisted | - Catgut Chromic |
| 8. SanaGut Plain | - Twisted | - Catgut Plain |
| 9. SanaLene | - Monofilament | - Polypropylene |
| 10. SanaLon | - Monofilament | - Polyamide |
| 11. SanaSil | - Black Braided | - Silk |
| 12. SanaBond | - Braided and Coated | - Polyester |
| 13. SanaWax | - Bone Wax 2.5g | |
| 14. SanaLene Mesh | - Knitted Polypropylene Mesh | |

Made in INDIA by:
Axa Parenterals Ltd

AXA House, Puhana Chowk, Roorkee, National Highway 73, Uttarakhand 247667
Please Visit us at www.sanacura.in | email : info@sanacura.in | Ph. : +91 72170 44614